

Intended use: Hydrozid® is used as cryotherapy for warts and other skin lesions.

Content: 200 ml HFC-134a (norflurane). Application Templates. Instructions for Use.

Instructions for Use: Hydrozid is applied to a lesion directly from the canister. Hydrozid should be applied from a distance of 2 to 3 cm and for 2 to 6 seconds, depending upon the location of the lesion. Hydrozid should be used according to illustration A, B, and C. In illustration A the spray head is locked. In illustration B the spray head is open and ready for use as in illustration C. Several repeated freeze-thaw cycles – up to 6 with 30 second intervals – are recommended as in a single treatment. The total duration of application to a lesion must be kept within the recommended maximum total (see panel). To avoid damaging the surrounding tissue use one of the application templates provided, carefully selecting the hole matching the size of the individual lesion. The spray jet intensity should also be adjusted to the selected hole in the application template. Some conditions may need repeated treatment to successfully remove a skin lesion. When repeat treatment is needed intervals of 1 to 2 weeks between treatments are recommended. Cryotherapy may produce a stinging or burning sensation on the skin. Necrosis of epidermal cells may result in the development of a blister. Full recovery takes about 10 to 14 days. Overexposure may lead to hypopigmentation, scarring or post-inflammatory hyperpigmentation.

Contraindications: Cryotherapy is contraindicated in patients with cryoglobulinaemia, Raynaud's disease, cold urticaria, blood dyscrasias, or uncertain diagnosis.

Caution: Never spray on mucous membranes, open wounds, cuts and delicate or eczematous skin. Caution is advised in terms of duration of treatment in children and people with sensitive and aging skin. Do not spray on healthy skin. Improper use may lead to unwanted damage to the skin and underlying tissues, including nerves. Do not inhale. Avoid contact with eyes. Keep out of reach of children. If medical advice is needed, have product container or label at hand.

Important: Hydrozid should only be used by healthcare professionals.

The application templates provided must be used only with Hydrozid. Due to the risk of infection, each application template should be used for one skin lesion only. Hydrozid is non-flammable and can be discarded in normal household waste.

Warning: Pressurised container; may burst if heated. Protect from sunlight. Do not expose to temperatures exceeding over 50 °C. Keep away from the heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Do not pierce or burn, even after use. Contains fluorinated greenhouse gas HFC-134a (0.248 kg 134a; 0.355 t CO₂ eq; GWP 1430).

Recommended maximum total treatment times for various types of skin lesion:

Skin Lesion	Maximum total treatment time in seconds
Verruca plantaris	36 (6x6seconds)
Verruca vulgaris	36 (6x6seconds)
Verruca plana	18 (6x3seconds)
Keratosis seborrhoica	36 (6x6seconds)
Lentigo (facial)	12 (6x2seconds)
Lentigo (non-facial)	36 (6x6seconds)
Acrochordon	36 (6x6seconds)
Keratosis actinica (facial)	12 (6x2seconds)
Keratosis actinica (non-facial)	36 (6x6seconds)
Condyloma	36 (6x6seconds)
Molluscum contagiosum	18 (6x3seconds)
Gingival melanin hyperpigmentation	36 (6x6seconds)
Cervical contact bleeding	36 (6x6seconds)

